

ST-120S Safety Enclosure

Mesa's newest feature to enhance operator safety



The ST-120S Open Top Safety Enclosure is the ideal solution for customers needing enhanced safety measures. Operators are protected from:

- Exposure to materials
- Shattering Glass
- Hazardous Projectiles
- Electrical and Mechanical Malfunctions

The Safety Enclosure is a robust 1" frame with ¼" thick clear abrasion resistant polycarbonate to protect the operator. The enclosure is also outfitted with magnetic proximity sensors that ensure operator compliance. The sensors are located on the frame as well as on the enclosure door to safeguard against operating the machine with the door open. Also, the ST-120S Safety Enclosure is equipped with a NEMA 4X rated push-button start switch that allows the operator to stand up to 5 feet away from the ST-120S cap torque analyzer.



For further information on the design and specifications of the safety enclosure, please contact torque@mesalabs.com or dial (303) 987-8000 to speak with a representative. We look forward to continuing to serve the cap testing industry and our loyal customer base.